BC 80

(D) Blutdruckmessgerät Gebrauchsanweisung	2-13
(GB) Blood pressure monitor Instructions for use	14-24
(F) Tensiomètre Mode d'emploi	25-36
(E) Tensiómetro Instrucciones de uso	37-48
(I) Misuratore di pressione Istruzioni per l'uso	49-60
(TR) Tansiyon ölçme cihazı Kullanım kılavuzu	61-71
(RUS) Прибор для измерения артериального давления на запястье Инструкция по применению	72-85
(PL) Ciśnieniomierz Instrukcja obsługi	86-97
Electromagnetic Compatibility Information	98-100

Inhalt

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Contents

1. Getting to know your device	14	6. Displaying and deleting measurements	22
2. Important notes	14	7. Transferring measurements	23
3. Unit description	17	8. Error messages/troubleshooting	23
4. Preparing for the measurement.....	18	9. Cleaning and maintenance.....	23
5. Measuring blood pressure.....	19	10. Technical specifications	24

Dear customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, gentle therapy, blood pressure/diagnosis, weight, massage, beauty and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the notes they contain.

With kind regards,
Your Beurer team

1. Getting to know your device

Check that the packaging of the Beurer BC 80 blood pressure monitor has not been tampered with and make sure that all the required contents are present.

The wrist blood pressure monitor is used to carry out non-invasive measurement and monitoring of arterial blood pressure values in human adults. This allows you to quickly and easily measure your blood pressure, save the measurements and display the development of the measurements. You are also warned of possible existing cardiac arrhythmia.

The recorded values are classified according to the WHO guidelines and evaluated graphically.


Store these instructions for use for future reference and make them accessible to other users.

2. Important notes

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF
	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Manufacturer
Storage 	Permissible storage temperature and humidity
Operating 	Permissible operating temperature and humidity
	Protect from moisture
SN	Serial number

The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are no substitute for a medical examination. Discuss the measurements with your doctor and never base any medical decisions on them (e.g. medicines and their administration).
- Do not use the blood pressure monitor on newborns or patients with pre-eclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- In the case of restricted circulation on the arm as a result of chronic or acute vascular diseases (including vascular constriction), the accuracy of the wrist measurement is limited. In this case you should avoid using an upper arm blood pressure monitor.

- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people who have the specified wrist measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- The blood pressure monitor can only be operated with batteries.
- To conserve the batteries, the monitor switches off automatically if you do not press any buttons for 2 minutes.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Instructions for storage and maintenance

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measurements and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
- Do not press the buttons before the cuff is placed on the arm.
- We recommend that the batteries be removed if the device will not be used for a prolonged period of time.

Notes on batteries

- Swallowing batteries can be extremely dangerous. Keep the batteries and products out of the reach of small children. If a battery is swallowed, seek medical attention immediately.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- If the batteries are empty or you do not intend to use the device for longer periods, remove the batteries from the device. This prevents damage which could occur when the batteries expire. Always replace all batteries at the same time.
- Do not use different types of batteries, different brands or batteries with different capacities. Use alkaline batteries.

i Instructions for repairs and disposal


- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not open the device. Failure to comply will result in voiding of the warranty.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers. Before making a claim, please check the batteries first and replace them if necessary.
- Please dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.


3. Unit description

1. WHO indicator
2. USB interface
3. Display
4. Positioning display (heart level)
5. Memory button **M1**
6. Memory button **M2**
7. **START/STOP** button **Ⓢ**
8. Battery compartment lid
9. Wrist cuff

Information on the display:


1. WHO classification
2. Time and date
3. Systolic pressure
4. Diastolic pressure
5. Resting indicator (symbol is inactive)
6. Cardiac arrhythmia symbol  Pulse symbol 
7. Calculated pulse value
8. User memory  / 
9. Pump up, release air (arrow)
10. Memory space number/memory display average value (A), morning (AM), evening (PM)
11. Battery status indicator 

4. Preparing for the measurement

Inserting the batteries

- Remove the battery compartment lid on the left side of the device.
- Insert two 1.5 V micro (alkaline type LR03) batteries. Make sure that the batteries are inserted the correct way round. Do not use rechargeable batteries.
- Close the battery compartment lid again carefully.


If the battery change symbol  is flashing and **E6** appears, no more measurements are possible and you must replace all batteries.

All display elements are briefly displayed, **24H** flashes in the display. Now set the date and time as described below.

Used batteries should not be disposed of in normal household waste. Dispose of them via your electronics retailer or your local recycling point. You are legally obligated to do so.

Note: The codes below are printed on batteries containing harmful substances: Pb: battery contains lead, Cd: battery contains cadmium, Hg: battery contains mercury.

Setting the hour format, date and time

This menu allows you to set the following functions, one after another.

Hour format

Date

Time

It is essential to set the date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

❗ If you press and hold the **M1** or **M2** memory button, you can set the values more quickly.

Hour format

- Press and hold the **START/STOP** button ❶ for 5 seconds.
- Select the desired hour format using the **M1/M2** memory buttons and confirm with the **START/STOP** button ❶.

Date

The year flashes on the display.

- Select the year using the **M1/M2** memory buttons and confirm with the **START/STOP** button ❶.

The month flashes on the display.

- Select the month using the **M1/M2** memory buttons and confirm with the **START/STOP** button ❶.

The day flashes on the display.

- Select the day using the **M1/M2** memory buttons and confirm with the **START/STOP** button ❶.

❗ If the hour format is set as 12h, the day/month display sequence is reversed.

Time

The hours flash on the display.

- Select the hour using the **M1/M2** memory buttons and confirm with the **START/STOP** button ❶.

The minutes flash on the display.


- Select the minute with the **M1/M2** memory buttons and confirm with the **START/STOP** button ❶.

Once all the data has been set, the device switches off automatically.

5. Measuring blood pressure

Attaching the cuff

- Expose your left wrist. Ensure that the circulation of the arm is not hindered by tight clothing or similar. Place the cuff on the inside of your wrist.

- Close the cuff with the hook-and-loop fastener such that the upper edge of the device is positioned approx. 1 cm underneath the heel of the hand.
- The cuff must be fitted tightly around the wrist, but must not constrict.

Adopt the correct posture

- Before every measurement, relax for about five minutes. Otherwise deviations can occur.
- You can take the measurement while sitting or lying. To take your blood pressure, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground. Make sure to rest your arm and move it. Always make sure that the cuff is at heart level. Otherwise significant deviations can occur. Relax your arm and the palm of your hand. An O.K. symbol bar is integrated into the device to additionally assist you. It signals the correct position of the blood pressure monitor. You are in the right position if you can see O.K. on the display.
- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.


Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- Press the **START/STOP** button **ⓘ** to start the blood pressure monitor. All displays will illuminate briefly.

The blood pressure monitor will begin the measurement automatically after 3 seconds.


The measurement is taken during the inflation phase.

- You can cancel the measurement at any time by pressing the **START/STOP** button **ⓘ**.

As soon as a pulse is found, the pulse symbol **♥** is displayed.

- The systolic pressure, diastolic pressure and pulse rate measurements are displayed.

- **E_** appears if the measurement could not be performed properly. Observe the chapter on error messages/troubleshooting in these instructions for use and repeat the measurement.
- Now select the desired user memory by pressing the **M1** or **M2** memory buttons. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant symbol **M1** or **M2** appears in the display.
- Press the **START/STOP** button **ⓘ** to switch off the blood pressure monitor. The measurement is then stored in the selected user memory.

If you forget to turn off the device, it will switch off automatically after approx. 3 minutes. In this case too, the value is stored in the selected or most recently used user memory.

- Wait at least 5 minutes before taking another measurement.

Evaluating results

Cardiac arrhythmia:

This device can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol .

This can be an indicator for arrhythmia. Arrhythmia is an illness in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heartbeat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess stimulants, stress or lack of sleep. Arrhythmia can only be determined through an examination by your doctor.

If the symbol  is shown on the display after the measurement has been taken, repeat the measurement. Please ensure that you rest for 5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor. Self-diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

WHO classification:

In accordance with the guidelines/definitions of the World Health Organization and the latest findings, the measurements can be classified and assessed according to the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The bar chart on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different WHO categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical WHO classification on the device always shows the higher category; for the example given this would be 'High normal'.

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
Grade 3: severe hypertension	≥180	≥110	seek medical attention
Grade 2: moderate hypertension	160-179	100-109	seek medical attention
Grade 1: mild hypertension	140-159	90-99	regular monitoring by doctor
High normal	130-139	85-89	regular monitoring by doctor
Normal	120-129	80-84	self-monitoring
Optimal	<120	<80	self-monitoring


Source: WHO, 1999

6. Displaying and deleting measurements

User memory

The results of every successful measurement are stored together with the date and time. If there are more than 60 measurements, the oldest measurements are lost.

- To access memory recall mode, press the **M1** or **M2** memory button to select your desired user memory (1 2).
- To view the measurements for user memory 1, press the **M1** memory button.
- To view the measurements for user memory 2, press the **M2** memory button.

Your last measurement will appear on the display.


Average values

- Press the relevant memory button (**M1** or **M2**).
- ① Press the **M1** button if you have selected user memory 1. Press the **M2** memory button if you have selected user memory 2.

A flashes on the display.

The average value of all saved measurements in this user memory is displayed.


- Press the relevant memory button (**M1** or **M2**).

A flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).

- Press the relevant memory button (**M1** or **M2**).


P flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

Individual measured values

- When the relevant memory button (**M1** or **M2**) is pressed again, the last individual measurement is displayed (in this example, measurement 03).
- When the relevant memory button (**M1** or **M2**) is pressed again, you can view your individual measurements.
- To switch the device off again, press the **START/STOP** button ①.

- ① You can exit the menu at any time by pressing the **START/STOP** button ①.

Deleting measured values

- To clear the memory of the relevant user memory, you must first select a user memory.

- Start individual measurement access.
- Press and hold the **M1** & **M2** memory buttons for 5 seconds.

CL and **00** appear on the display.

All the values in the current user memory are deleted.

The device will then switch off automatically.

7. Transferring measurements

Connect the blood pressure monitor to your PC using the USB cable.

- ⓘ No data transfer may be launched whilst performing a measurement.

PC is shown on the display. Begin the data transfer in the “HealthManager” PC software. During the data transfer, an animation is shown on the display. A successful data transfer is displayed as in figure 1. If the data transfer is unsuccessful, an error message appears as in figure 2. In this case, interrupt the PC connection and start the data transfer again.


figure 1

figure 2

After 30 seconds of not being in use or if communication with the PC is interrupted, the blood pressure monitor switches itself off automatically.

8. Error messages/troubleshooting

In the event of errors, the error message **E_** appears on the display.

Error messages may appear if:

- It was not possible to record the pulse: **E1**;
- You move or speak during the measurement: **E2**;
- The cuff is fastened too tightly or loosely: **E3**;
- Errors occur during the measurement: **E4**;
- The pump pressure is higher than 300 mmHg: **E5**;
- The batteries are almost empty  **E6**.

In such cases, repeat the measurement. Ensure that you do not move or speak. If necessary, reinsert or replace the batteries.

9. Cleaning and maintenance

- Clean your blood pressure monitor carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances should you hold the unit under water, as this can cause liquid to enter and damage the unit.
- Never place any heavy objects on the device.

10. Technical specifications

Model no.	BC 80
Measurement method	Oscillometric, non-invasive blood pressure measurement on the wrist
Measurement range	Cuff pressure 0-300 mmHg, Systolic 50-250 mmHg, Diastolic 30-200 mmHg, Pulse 40-180 beats/minute
Display accuracy	Systolic ± 3 mmHg, Diastolic ± 3 mmHg, Pulse ± 5 % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: Systolic 8 mmHg / Diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 97 mm x W 68 mm x H 22 mm
Weight	Approximately 108 g (without batteries)
Cuff size	135 to 230 mm
Permissible operating conditions	+5 °C to +40 °C, 15-93 % relative humidity (non-condensing)
Permissible storage conditions	-25 °C to +70 °C, 10-93 % relative humidity, 700-1050 hPa ambient pressure
Power supply	2 x 1.5V --- AAA batteries
Battery life	For approx. 300 measurements, depending on levels of blood pressure and pump pressure

Accessories	Instructions for use, 2 x 1.5V AAA batteries, storage box
Classification	Internal supply, IPX0, no AP or APG, continuous operation, application part type BF

- This device is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Services address or found at the end of the instructions for use.
- This device corresponds to the EU Medical Devices Directive 93/42/EC, the German Medical Devices Act (Medizinproduktegesetz) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

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Electromagnetic Compatibility Information

Guidance and manufacture's declaration – electromagnetic emissions – for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission

The BC80 blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the BC80 blood pressure monitor should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The BC80 blood pressure monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The BC80 blood pressure monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacture's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity

The BC80 blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the BC80 blood pressure monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacture's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacture's declaration – electromagnetic immunity

The BC80 blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the BC80 blood pressure monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BC80 blood pressure monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of  equipment marked with the following symbol:</p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BC80 blood pressure monitor is used exceeds the applicable RF compliance level above, the BC80 blood pressure monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the BC80 blood pressure monitor.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for

EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING*Recommended separation distances between portable and mobile RF communications equipment and the BC80 blood pressure monitor*

The BC80 blood pressure monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BC80 blood pressure monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BC80 blood pressure monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.117	0.117	0.234
0.1	0.370	0.370	0.740
1	1.170	1.170	2.340
10	3.700	3.700	7.400
100	11.7	11.7	23.4

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.